

# IBUPROFEN

# Ibuprofen trial set-up

## Pre-blend (before PD-Granulation)

- 95.0% Ibuprofen
- 5.0% Ac-Di-Sol


## Final blend\*

- 95.7 % Pre-blend
- 4.0 % Klucel EXF
- 0.3 % Aerosil 200

Final tablet contains 90.9 % of Ibuprofen

\* Externally lubricated with Magnesiumstearate

# COMPACTIBILITY of PDG-Ibuprofen



# Properties of PDG-Ibuprofen tablets

aspect	
Crushing Force (9mm concave tablets (rc=15 mm), 312.5 mg)	104 N
Tablet Strength	1.66 MPa
Friability	0,32 %
Disintegration	< 15 min
Dissolution	not yet tested, but is not expected to be a problem

simulated tableting speed: 43000 tablets/hr (Fette 102i, 30 rpm), no sticking to punches observed

# Conclusion PDG Ibuprofen

So, a rather difficult-to-process material like Ibuprofen results in

**a well flowable granulate/tableting mix**

and tablets, which

**are strong enough:** 104N (= 1.66 MPa), low friability, no capping, no sticking to punches

**disintegrate within:** 15 minutes

**dissolution:** probably OK

**have high drug load:** 90.9 %

## Pneumatic Dry Granulation results in

### Avicel PH105

- ❖ excellently flowable granules
- ❖ obtained by roller compaction to low density ribbons
- ❖ thus ensuring optimal recompaction properties

### Paracetamol Ibuprofen

- ❖ enables tablet formulations with difficult to compact API's
  - hardness OK, friability OK, no capping
  - disintegration and dissolution OK
  - high drug loads

Virtually every API, that in principle can be compacted, can also be  
**PD-G**ranulated !

The result is a flowable and compactable granulate !

This extends the application of dry granulation enormously:

**PDG**-technology now provides you with  
an excellent alternative to wet granulation  
without the disadvantages and  
being far more economical

**Contact us at [www.atacamalabs.com](http://www.atacamalabs.com)  
for more results**